

TECHNICAL DATA SHEET

TEX Q®

58 Kill Claims – One Step Cleaner and Disinfectant

TECHNICAL DATA SHEET

Quaternary Ammonium Disinfectant

TexQ® Disinfectant is the last generation of quaternary ammonium compounds (QACs), EPA registered one-step cleaner/disinfectant. It is effective against a broad spectrum of bacteria, viruses and fungi, and inhibits the growth of mold and mildew and their odors when used as directed.

TexQ® Disinfectant is available as Ready-to-Use solution in a spray bottle (reorder number TX650) and Ready-to-Use solution in one gallon bottle (reorder number TX652). All solutions are 0.2 micron filtered.

Industries

- Biotechnology
- Hospitals, Pharmacies
- Medical device manufacturing facilities
- Pharmaceutical manufacturing facilities
- Veterinary clinics and laboratories
- Food processing facilities

Applications

- Cleaning and disinfecting hard, non-porous surfaces
- Cleaning and disinfecting small surfaces (tables, equipment, isolators, hoods)
- Cleaning and disinfecting large surfaces (floors, walls, ceilings)
- Recommended for use as part of a disinfectant rotation program

Surfaces

- Plastic
- Polycarbonate
- Plexiglas®
- Glass
- Chrome
- Vinyl Upholstery
- Stainless Steel



Made in USA

Products

Number	Description	Packaging
TX650	TexQ® Disinfectant Ready-to-Use (RTU)	22 oz. trigger spray, 12 bottles/case, gamma-irradiated
TX652	TexQ® Disinfectant Ready-to-Use (RTU)	One gallon concentrate, 4 bottles/case, gamma-irradiated



Features & Benefits

- **58 Kill Claims** – Proven efficacy against the most common cleanroom bacteria, viruses, and fungi
- **The Ready-to-Use solutions are Gamma-irradiated** – Compliance with aseptic environment requirements and USP <797>
- **Complex Formulation** – Excellent cleaning and disinfecting properties for uniform and complete disinfection of surface
- **One-step cleaner/disinfectant** – Cleans and disinfects in one step. Saves additional costs for the cleaner. Saves time and labor. Easy to use
- **Free of dyes and fragrance** – No respiratory irritation – safe for the staff, no additional contamination from the vapors
- **Functional use label on bottle** – Easy documentation and usage control, record the date opened and operator initials
- **EPA registered** – Kill claims provide assurance of efficacy of disinfection

TECHNICAL DATA SHEET

TexQ® Disinfectant Kill Claims

TX650 TexQ® Disinfectant	Contact Time in Minutes*	TX650 TexQ® Disinfectant	Contact Time in Minutes*
Bacteria		Viruses	
Community Associated Methicillin Resistant Staphylococcus Aureus	10	Avian influenza A Virus (<i>H5N1</i>)	10
Methicillin Resistant Staphylococcus Aureus	10	Avian Influenza A/Turkey/Wisconsin	10
Burkholderia cepacia	10	Canine Coronavirus	10
Campylobacter jejuni	10	Canine Distemper Virus	10
Corynebacterium ammoniagenes	10	Duck Hepatitis B Virus	10
Enterobacter aerogenes	10	Hantavirus	10
Enterobacter cloacae	10	Hepatitis B virus (<i>HBV</i>)	10
Enterobacteriaceae w/extended beta lactamase resistance	10	Hepatitis C virus (<i>HCV</i>)	10
Enterococcus faecalis	10	Herpes Simplex Type 1	10
Enterococcus faecium (<i>Vancomycin resistant</i>)	10	Herpes Simplex Type 2	10
Escherichia coli	10	HIV-1 (<i>AIDS virus</i>)	2
Escherichia coli (<i>Antibiotic resistant</i>)	10	Human Coronavirus	10
Escherichia coli O157:H7	10	Infectious Bovine Rhinotracheitis virus (<i>IBR</i>)	10
Klebsiella pneumoniae	10	Influenza Type A / Brazil	10
Klebsiella pneumoniae (<i>Antibiotic resistant</i>)	10	Influenza A (H1N1) Virus	10
Legionella pneumophila	10	Newcastle Disease virus	10
Listeria monocytogenes	10	Porcine Respiratory & Reproductive Syndrome Virus (<i>PRRSV</i>)	10
Proteus mirabilis	10	Porcine Rotavirus	10
Proteus vulgaris	10	Pseudorabies Virus (<i>Rabies Virus</i>)	10
Pseudomonas aeruginosa	10	Respiratory Syncytial Virus (<i>RSV</i>)	10
Pseudomonas aeruginosa (<i>Antibiotic resistant</i>)	10	Transmissible Gastroenteritis (<i>TGE</i>)	10
Salmonella enterica	10	Vaccinia Virus	10
Salmonella schottmuelleri	10		
Salmonella typhi	10	Fungi	
Serratia marcescens	10	Candida albicans	10
Shigella dysenteriae	10	Trichophyton mentagrophytes (<i>Athlete's Foot Fungus</i>)	10
Shigella flexneri	10		
Shigella sonnei	10	Mildewstatic	
Staphylococcus aureus	10	Aspergillus niger	
Staphylococcus epidermidis (<i>Antibiotic resistant</i>)	10		
Streptococcus pyogenes	10		
Vibrio cholerae	10		

*As of 10/2017. Subject to change per EPA review.

†Tested according to the AOAC Use Dilution test method on hard inanimate surfaces, in the presence of 5% organic serum.

Shelf life:

For TX650 TexQ® Disinfectant Ready-to-Use (RTU), both 22 oz and 1-gallon bottles:

2 years (the expiration date is indicated on the product label), even after opening the bottle or spraying from the bottle.

Disposal:

Dispose of contents and container in accordance with all local, regional, national and international regulations (see Safety Data Sheet, Section 13.)



Americas	<i>+1 336 996 7046</i>
EMEA	<i>+31 88 1307 410</i>
Philippines	<i>+63 49 543 0241</i>
China	<i>+86 512 6303 3700</i>

texwipe.com